



The EarPopper® middle ear inflation device balances pressure in the middle ear by delivering a safe, constant stream of air into the nasal cavity. Treats common middle ear pressure problems without medication or surgery like:

- Eustachian tube dysfunction
- Otitis media with effusion (glue ear)
- Ear fullness caused by air travel or scuba diving

The EarPopper's constant stream of air opens the Eustachian Tube during the moment of swallowing, clearing the middle ear, relieving negative ear pressure and allowing any fluids to drain. The EarPopper is suitable for adults and children 4+ (with adult supervision).

*USA Ordering Information (prescription required)*

PRODUCT NO.	DESCRIPTION	UNITS/BOX
EP-2100	Includes 4 AAA Batteries	1
EP-1010	Disposable Tips	50

*Rest of World Ordering Information*

PRODUCT NO.	DESCRIPTION	UNITS/BOX
EP-2100A	Includes 4 AAA Batteries	1

*Clinically proven to relieve negative ear pressure.*



Flying



Fluid in Ear



Diving & Swimming

Consult a physician and read the instructions found in package before using the EarPopper. Do not use the EarPopper if you have an ear infection, upper respiratory infection, nasal congestion, cold symptoms or a perforated eardrum. The EarPopper is not intended for multiple users.



POINT CAMERA TO READ IFU

M-77 | RevA | OCT2021